FR

FORMULE

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PROPERTIES

PHARMAETHYL has a greater cryo-anaesthetic activity than ethyl chloride thanks to its lower boiling point. When PHARMAETHYL is sprayed as a fog made of suspended microdrops, it evaporates faster than ethyl chloride and **its cooling action is quicker and deeper.**

INDICATIONS

Non-flammable cryo-anaesthetic, to be used in local application on the mucosa of the oral cavity and its adnexa: lancing of abscesses, extraction of deciduous teeth or pyorrheal tooth roots, vitality test.

INSTRUCTIONS FOR USE

Remove the white cap from the top of the bottle. Take the canula out of the cellophane bag and fit it on the end of the tube by pressing it home.

To spray, exert **vertical** pressure on the pushbutton strong enough to deliver one dose of product. Spray a 2 to 4 cm distance away from the point to be anaesthetised.

To obtain a wider cryo-anaesthesia, deliver the required doses, taking care not to target always the same point, to prevent any risk of burns. After each use, remove the canula, clean it and disinfect it in a suitable disinfecting solution.

SPECIAL WARNINGS AND PRECAUTIONS FOR USE

Pressurized container: protect from sunlight and do not expose to temperatures exceeding 50°C. Do not pierce or burn, even after use.

Do not spray on a naked flame or any incandescent material.

Keep away from sources of ignition - No smoking.

Xi – Irritant (contain (S)-p-mentha-1,8-diene)

May cause sensitisation by skin contact.

Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

Do not breathe vapour.

Avoid contact with skin. Wear suitable gloves.

Use only in well-ventilated areas.

Avoid release to the environment. Refer to special instructions/safety data sheets.

PHARMAETHYL is not toxic however, it is **designed for local anaesthesia only.**

PRESENTATION

1 x 150 ml aerosol can + 1 canula.

Clinical presentation: 3 aerosol cans.

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